




TIME MEDICAL
SYSTEMS



PANION
SUPREME
veterinary MRI system



PANION SUPREME Features



Dry magnet - forget about helium, PANION SUPREME is a cryogen-free MR.



No quench pipe required - significantly reducing construction costs of new installations.



Strongest Gradients - 52 mT/m with slew rate at 142 mT/m for ultra-high quality scans.




PRODIVA drives user-friendly clinical workflows during patient set-up, scanning and data management.




Widest Range of Clinical Applications including: ce-MRA, DIXON Fat & Water Separation, DWI & ADC map.


Clinical Images




Brain T2 FSE Sagittal



C-Spine T2 FSE Sagittal




Brain T2 FSE Dorsal




Brain T1 FSE Sagittal Post Contrast


Clinical Images




Brain + C-Spine T2 FSE Sagittal



C-Spine T2 FSE Axial




TL-Spine IRFSE Sagittal




TL-Spine T2 FSE Dorsal


Clinical Images




TL-Spine T2 FSE Axial



TL-Spine T2 FSE Axial



Renal T2 FSE Dorsal





Crab T2 FSE Dorsal
(1024 matrix high resolution)



Remote Support Services

- Provide real-time remote software support through internet from authorized supporting personnel.
- Provide real-time monitoring of system hardware status from anywhere.
- Provides insight for preventive maintenance.
- Efficient service turnaround time when service engineers know the hardware status from the service centres remotely.



Prodiva Cloud Teleradiology Platform

- Images sent to the Prodiva Cloud for secure storage and review.
- Full-featured web-based interface allows diagnosis and reporting by remote experts anywhere, anytime without boundary.
- Intelligent deep learning engine provides first-pass diagnostic assistance immediately upon arrival of images (WIP)



Time Medical USA

890 Cowan Road, Suite J, Burlingame,
CA 94010, USA

Time Medical Hong Kong

Rm. 301, Building 20E, Science Park E Ave.
Hong Kong Science & Technology Park,
Shatin, N.T. Hong Kong

Time Medical Singapore

Blk 71, Ayer Rajah Crescent
#03-04/05
Singapore 139951

Time Medical Shanghai

4/F, Block B, 563 Song Tao Road,
Zhangjiang Hi-Tech Park,
Pudong District,
Shanghai, 201203,
China

Time Medical Jiangsu

G18 China Medical City,
Taizhou,
Jiangsu Province, 225300,
China

sales@time-medical.com | www.time-medical.com

Until Time Medical releases these products for general sale, Time Medical reserves the right to alter prices, specifications, features, capabilities, functions, manufacturing release dates, general availability of the products at any time.